



Sani-Cloth® AF³

Cleaning & Disinfectant



Contents


Healthcare Associated Infections: A Threat to Patient Safety.....	4
Why Choose AF ³	6
Antimicrobial Efficacy.....	8
Product Information.....	9
Point of Use Accessories.....	9
How to Use.....	10
Iconography Explained.....	11
Wipe Selector Guide.....	12
Training & Educational Resources.....	13
Made in Britain.....	13



Healthcare Associated Infections

A Threat to Patient Safety

Healthcare acquired infections (HAIs) in hospitals alone cause more deaths in Europe than any other infectious disease under surveillance at the European Centre for Disease and Control (ECDC).



Common HAIs

Respiratory Tract Infections (RTI)
Example: *Pneumonia*

Surgical Site Infections (SSI)
Example: *Superficial and deep SSI*

Urinary Tract Infections (UTI)
Example: *Catheter-associated urinary tract infection (CAUTI)*

Bloodstream Infections (BSI)
Example: *Escherichia coli (E.coli)*

Gastrointestinal Infections (GI)
Example: *Clostridium difficile*

Checklist to Prevent HAIs⁸

- ✓ Infection prevention and control, including hand hygiene.
- ✓ Antimicrobial stewardship programmes.
- ✓ Healthcare worker training.
- ✓ Improved information for patients and visitors about how they can prevent HAIs and use antibiotics prudently.
- ✓ HAI surveillance at local and national level.
- ✓ Improved microbiological laboratory support.

Why Choose Sani Cloth® AF³?

Sani-Cloth® AF³ is an alcohol free, dual quaternary ammonium formulation, ideal to use for routine cleaning and disinfecting.

Effective Kill at Speed

Sani-Cloth® AF³ is clinically proven to be effective against bacteria, enveloped and non-enveloped viruses (including norovirus), and yeast in 60 seconds.




Free From

AF³ wipes are free from PHMB and pose no known risk of triggering allergic reactions.

Compatible by Design

AF³ has undergone rigorous ASTM and ISO material testing and has been approved with more than 450 medical devices.



Compliance Simplified

Pre-dosed, 2-in-1 disposable cleaning and disinfection wipes that require no mixing or measuring, making surface decontamination reliable, easy and efficient every time.

Made in Britain

Sani-Cloth® AF³ is manufactured in the UK, resulting in supply chain reliability and carbon footprint reduction.



Tried and Trusted

The number one alcohol free wipe in North America!



TUBERCULOCIDAL, VIRUCIDAL, YEASTICIDAL & BACTERICIDAL IN 1 MINUTE

ACTIVE INGREDIENTS:	
N-alkyl dimethyl benzyl ammonium chloride	0.14%
N-alkyl dimethyl ethylbenzyl ammonium chloride	0.14%
OTHER INGREDIENTS	99.72%
TOTAL	100.00%

Does not include the weight of the wipe

CONTAINS: 200 Wipes
200mm x 200mm



Antimicrobial Efficacy

Sani-Cloth® AF³ has been independently tested under rigorous EN standards. Testing is performed using liquid extracted from the wipe under dirty conditions, representing real-life situations more accurately.

Fully bactericidal, virucidal, yeasticidal, and tuberculocidal in 1 minute.

Test Method		Contact Time	Conditions
Bacteria			
EN 16615	Pseudomonas aeruginosa	1 min	Dirty
EN 16615	Staphylococcus aureus	1 min	Dirty
EN 16615	Enterococcus hirae	1 min	Dirty
EN 13727	Pseudomonas aeruginosa	1 min	Dirty
EN 13727+B8:B15	Staphylococcus aureus	1 min	Dirty
EN 13727	Enterococcus hirae	1 min	Dirty
EN 13727	MRSA	1 min	Dirty
EN 13727	E.coli	1 min	Dirty
Yeast			
EN 16615	Candida albicans	1 min	Dirty
EN 13624	Candida albicans	1 min	Dirty
EN 13624	Candida auris	1 min	Dirty
Enveloped Virus			
EN 14476	Vaccinia virus	1 min	Dirty
EN 14476	Mers-Cov	1 min	Dirty
EN 14476	HIV	1 min	Dirty
EN 14476	Hepatitis B	1 min	Dirty
EN 14476	Hepatitis C12:C19s C	1 min	Dirty
EN 14476	SARS-CoV-2 / Coronavirus	1 min	Dirty
EN 14476	Influenza H1N1	1 min	Dirty
EN 14476	Influenza H3N2	1 min	Dirty
EN 14476	Influenza H5N1	1 min	Dirty
EN 14476	Measles	1 min	Dirty
EN 14476	Mumps	1 min	Dirty
Non-Enveloped Virus			
EN 14476	Norovirus	1 min	Dirty
TB			
EN 14348	Mycobacterium terrae	1 min	Dirty

(c) clean conditions / (d) dirty conditions

For full testing data, please contact your local representative

Product Information

Sani-Cloth®
CLEANING AND DISINFECTANT

AF³
LARGE WIPE

CE
2797

DISINFECTS IN
1
MINUTE

PHMB
FREE

COMPATIBLE
BY DESIGN

MADE IN
BRITAIN

Sani-Cloth® AF³ wipes are alcohol free, ready to use low linting, disposable 2 in 1 cleaning and disinfectant wipes for non-invasive medical devices and general surfaces. Bactericidal, virucidal, yeasticidal, and tuberculocidal in 60 seconds..

	WIPE SIZE	WIPE COUNT	QTY PER CASE	CLASSIFICATION	REORDER NO	NHS CODE	ROUTE
200ct Flow Pack	267mm x 210mm	200	6	MD Class IIa / Biocide	EAF3001	VJT592	NHS Supply Chain

Point of Use Accessories

Our point of use accessories help put Super Sani-Cloth® in the right place at the right time to help you do your job more efficiently and effectively.

	ACCESSORY	REORDER NO	NHS CODE	ROUTE
	Canister Bracket Our canister bracket affixes securely to the wall allowing easy, one-handed wipe dispensing.	A001GRY	N/A	Direct / Distributor
	Soft Pack Dispenser Our Soft Pack Dispenser has been designed to hold the lid of the flow wrap securely in place while in use, facilitating convenient single-handed dispensing of wipes. Additionally, the dispenser incorporates antimicrobial additives that have been validated to decrease bacteria growth by up to 99.99% (ISO 22196:2011).	A004GRY	N/A	Direct / Distributor
	I Am Clean Sticky Notes Our recyclable sticky notes can be used to easily and effectively identify equipment or surfaces that have been cleaned.	A001YEL	N/A	Direct / Distributor
	I Am Clean Sticky Tape Our tape can be dispensed to whatever length is required and is used to easily identify and communicate that a surface or piece of equipment has been cleaned.	A001RED	N/A	Direct / Distributor

How to Use

STEP 1: Risk Assessment



Follow your local protocol for personal protective equipment (PPE), hand hygiene and disposal guidelines before using the wipes.

STEP 2: Cleaning & Disinfecting



If the surface is visibly soiled, it is recommended to clean the surface prior to using a fresh wipe to disinfect.



Wipe in an “S-shaped” motion from top to bottom



Ensure the surface remains visibly wet for 1 minute for complete disinfection. Use an additional wipe if necessary. Allow the surface to air dry.

STEP 3: Disposal



Do not re-use wipes.



Do not flush. Discard used wipes in the appropriate waste bin following your local policy.

Recommended Uses



High touch point surfaces



Common equipment and large surface areas



Iconography Explained

We believe total transparency to be essential for earning the trust of our customers. We take pride in ensuring absolute clarity on product labels, enabling healthcare professionals to carry out their duties with increased efficiency and, above all, safely.

Icon	Description
	Product name and descriptor that this is a 2-in-1 cleaning and disinfectant wipe.
	Complete disinfection of the listed pathogens on the packaging in 1 minute.
	Does not contain PHMB.
	Do not flush.
ACTIVE INGREDIENTS	* n-Alkyl Dimethyl Benzyl Ammonium Chloride...0.14% * n-Alkyl Dimethyl Ethylbenzyl Ammonium Chloride...0.14%
	Period after opening – this is the period of time the product can be used after opening. (Not to be confused with an expiry date which relates to an unopened product).
	Indicates that the medical device (wipe) is intended for one use.
	Products packaged in various materials. Refer to packaging for recyclability details and check with local policy disposal guidelines.
	Indicates the medical device manufacture's address: PDI Ltd, Pywell Road, Willowbrook East Ind Estate, Corby, NN17 5XJ UK
	European Authorised Representative: Nex Medical Antiseptics Srl, Via Per Arluno 37/39-20010 Casorezzo (MI) Italy.
	The UKCA marking indicates conformity with the applicable requirements for products sold within Great Britain.
	CE mark indicating the product is a class IIa medical device.

Wipe Selector Guide

We understand that there are different variables and sensitivities in every healthcare setting. Therefore, as part of our layered approach to infection prevention, we assess customers current practices and challenges before proposing the solution that help to make their current practices easier and more efficient.

						
Performance						
Cleaning	✓	✓	✓		✓	✓
Disinfecting	✓	✓	✓	✓		
Use						
Routine between patients		✓	✓	✓	✓	✓
Outbreak prevention & management	✓	✓	✓	✓		
C.diff Outbreak	✓					
Formulation						
Sodium Hypochlorite	✓					
Dual quat + alcohol		✓				
Quat (alcohol-free)			✓			
Alcohol				✓		✓
Detergent					✓	
Efficacy						
Sporicidal - C.diff	✓					
Tuberculocidal		✓	✓	✓		
Virucidal - Norovirus	✓	✓	✓			
Virucidal - Enveloped Viruses	✓	✓	✓	✓		
Yeasticidal	✓	✓	✓	✓		
Bactericidal	✓	✓	✓	✓		
Contact time						
Total contact time of listed pathogens	2 minutes	1 minute	1 minute	1 minute	N/A	N/A

Training & Educational Resources

It doesn't matter how good a product is; if it's not used in the correct manner it will not be effective at doing its job and could result in putting patients, healthcare workers and visitors at risk of infection.

Therefore, as part of our layered approach to infection prevention, not only do we assess and put forward solutions to meet your particular needs, we have a plethora of training and educational resources available to ensure you use our products correctly. Every time.

Training Resources

- ✓ IFU poster
- ✓ IFU range poster
- ✓ IFU video
- ✓ SWITCH posters
- ✓ FAQs

Educational Resources

- ✓ Bitesize learning series
- ✓ CPD validated courses

Made in Britain



Thankyou for choosing Made in Britain

Our entire range of Sani-Cloth® cleaning and disinfectant wipes are manufactured at our production site in Corby, Northamptonshire and as a result we have been officially granted Made in Britain membership.

What does this mean for you?

- You are helping reduce your carbon footprint.
- Our products are synonymous with quality, ethical business practices and trust.
- Reliable supply chain with no unexpected surcharges.
- You are supporting our economy and local employment.



Visit our factory!

If you would like to see where your wipes are manufactured, contact us today to arrange a visit!



References

1. <https://www.ecdc.europa.eu/sites/default/files/images/healthcare-associated-infections-threat-patient-safety-infographic.png>
2. <https://www.ecdc.europa.eu/sites/default/files/images/healthcare-associated-infections-threat-patient-safety-infographic.png>
3. <https://www.ecdc.europa.eu/sites/default/files/images/healthcare-associated-infections-threat-patient-safety-infographic.png>
4. Stewart et al. (2021) Impact of healthcare-associated infection on length of stay Journal of Hospital Infection 114 (2021) 23-31
5. <https://www.ecdc.europa.eu/sites/default/files/images/healthcare-associated-infections-threat-patient-safety-infographic.png>
6. Guest et al. (2019) Modelling the costs and consequences of reducing healthcare-associated infections by improving hand hygiene in an average hospital in England. BMJ Open 9:029971. doi:10.1136/bmjopen-2019-029971
7. de Kraker ME, Stewardson AJ, Harbarth S. Will 10 Million People Die a Year due to Antimicrobial Resistance by 2050? PLoS Med. 2016 Nov 29;13(11):e1002184. doi: 10.1371/journal.pmed.1002184. PMID: 27898664; PMCID: PMC5127510.
8. <https://www.ecdc.europa.eu/sites/default/files/images/healthcare-associated-infections-threat-patient-safety-infographic.png>



PDI EMEA Ltd. Pywell Rd, Willowbrook East Industrial Estate, Corby NN17 5XJ UK



+44 (0) 8081 697 945



contact@pdi-emea.com



pdi-intl.com