Compact and portable interpretive 12 lead ECG machine with 5" colour display



seca CT8000i-2

- + Quick and easy to use ECG machine
- + Compatible with all leading EMR systems emis, SystmOne, Vision and others
- + On-screen preview of ECG review recorded ECG prior to printing, saves paper
- + Outstanding signal quality for adult and paediatric ECG's due to high sample frequency
- + Advanced clinical interpretation ETM algorithm adult & paediatric
- + Compact lightweight design: 1.9 kg
- + Sealed keyboard ensures infection control
- + 5" colour LCD with lead off indicator





SCM (o) h(e

Compatible with emis and all EMR systems

seca CT8000i-2

Compact and portable interpretive 12 lead ECG machine with 5" colour display





The quick and easy ECG solution

The seca CT8000i-2 is the 2nd generation of our best-selling 12 lead ECG machine with an enhanced signal quality and new interpretation algorithms. It is compact, portable and has a user-friendly design, that allows for quick, easy and reliable ECG acquisition. Just enter patient data, control signal quality, acquire ECG and print it with the integrated high-resolution thermal printer. Taking an ECG has never been easier.

12

1.9 kg

5" colour high resolution LCD display

Storage of up to 100 ECGs

285 x 189 x 61 mm

Mains adapter

+ Print function

+ Save function

Rechargeable batteries,

+ Built-in thermal printer

+ Preview screen feature + Optimal signal quality

+ QTc interval - Ideal for Mental Health

+ QTc formulas = Bazett, Fredericia, Fremingham or Hodges

+ Compatible with all EMR systems emis, SystmOne, Vision and others + Sealed keyboard ensures infection

The seca CT8000i-2 makes ECG-reading very easy and comfortable for medical experts. The real-time signal preview and quality check with lead reversal detection guarantees proper electrode placement. Accommodating your preference, you have multiple ways to process the ECG: Print it with the build-in printer on roll or z-folded paper. Transfer it into your EMR system, export it as an PDF or save it on the device. The software provides a great range of QTc formulas like Bazett, Fredericia, Fremingham or Hodges.



+ 802Z.CT8000i-2: Z-fold ECG paper for seca CT8000i-2
+ 700.8000i-2: Carry case for seca CT8000i-2
+ seca trace Q-tab electrodes, pack of 500
+ Height Adjustable ECG Cart seca CT500-2
+ seca CT8000i-2-EX-W: extended warranty, giving 5 years piece of mind
Patient cable, mains cable, 300x seca trace Q-tab electrodes, alligator clips, pack of z-fold paper

united kingdom		
Barn Street / Digbeth	Т	+44 121 643
ingham B5 5QB	F	+44 121 633
ed Kinadom	F	info.uk@seca

Display

Memory

ECG Channels

Dimensions (W x H x D)

Net weight

Power supply

Functions

3 9349 3 3403

control

info.uk@seca.com

201115205021