

## Accessory Specification Sheet

### Safe Clip Wall Bracket

Product Code:

**45000000**

The Sharpsafe Near Patient Sharps Disposal (NPSD) system is a range of bracket, trays and trolleys which allow for sharps to be taken directly to the patient and disposed of at the point of use.

The Sharpsafe Safe Clip Wall Bracket is compatible with a range of Sharpsafe containers, from 1L - 9L High. The bracket allows for convenient and easy access to containers, to help minimise the risk of sharps injuries by allowing for disposal at the point of use.

#### Product Details

Product Dimensions (mm):	40 (L) x 74 (W) x 155 (H)
Product Weight:	35g
Material:	Recycled polypropylene
Compatible with Sizes:	1L, 1.8L, 2L, 3L, 4L, 7L and 9L high
GTIN:	5060107831209
Medical Device:	No
Reusable	Yes
Quantity per Carton:	10
Pack Dimensions (mm):	300 (L) x 240 (W) x 75 (H)
Pack Weight:	409g
Contains PVC:	No
Contains Latex:	No



FOR MORE INFORMATION:

✉ [sharpsafesales@vernagroup.com](mailto:sharpsafesales@vernagroup.com)

🌐 [vernacare.com](http://vernacare.com) ☎ +44 (0)1495 235800

Vernacare Limited, Block G, Newbridge Road Industrial Estate, Blackwood, Caerphilly, NP12 2YN, UK  
© Copyright Verna Group 2023. All rights reserved. GP0146\_A29\_01